

#### General Research Tabletop Oxygen Analyzer Item # 07-0006 Model O2

The Model O2 is a versatile research grade oxygen analyzer for measuring 5-100% oxygen concentrations. The unit includes an integrated gas sampling pump, flow controller, bright LED display for oxygen concentration and status indication, touch panel controls for ease of operation, Includes XC/LN Modes for 1 Volt DC analog and a terminal strip with relay outputs and a 4-20 ma analog output, a separate digital RS-232 output, and LEUR type front panel connection for response times. The front panel allows control of flow, alarm set-points, calibration, and other user features. The Model O2 has a CE clearance and can be used in research and other applications, but not approved for human medical use or diagnostics. Oxigraf OxiSoft software is also included with the unit for computer display and data logging.





# **Features**

## Precision Oxygen Analysis for Research and Industrial Applications

The Model O2 is a research-grade, single-channel tabletop oxygen analyzer designed for high-precision measurement of oxygen concentrations ranging from 5% to 100%. With an integrated gas sampling pump and advanced laser diode technology, this unit provides fast, reliable, and accurate oxygen analysis for laboratory, industrial, and safety applications.

## XC and LN Modes for Enhanced Accuracy

The Model O2 includes XC (Cross Calibration) Mode and LN (Laser Normalization) Mode, providing flexibility for different measurement conditions:

- LN Mode (Laser Normalization Mode): Optimized for nitrogen-based gas mixtures, this mode provides highprecision readings with a stability of ±0.2%, ensuring minimal drift in controlled environments.
- XC Mode (Cross Calibration Mode): Designed for applications involving complex gas mixtures, including noble gases, hydrocarbons, fluorocarbons, CO<sub>2</sub>, and N<sub>2</sub>O, this mode ensures accurate oxygen readings even in varying background gas compositions, with a stability of ±0.4%.

#### Fast & Reliable Oxygen Measurement

The Model O2 utilizes laser diode absorption technology, delivering precise oxygen concentration readings with an accuracy of  $\pm 0.2\%$  (LN mode) and a rapid response time of 150 milliseconds at a 250 mL/min flow rate. Unlike electrochemical sensors, it provides consistent performance without frequent calibration or sensor replacement, reducing downtime and maintenance costs.

### Integrated Sampling Pump for Continuous Monitoring

A built-in sampling pump ensures constant gas flow and accurate readings, making it ideal for laboratory research, industrial gas analysis, and safety monitoring applications. The unit's microprocessor-controlled flow sensor maintains stability, while a hydrophobic filter prevents contamination and extends the device's lifespan.

#### **Multiple Output Options for Seamless Integration**

The Model O2 supports analog (0-1V DC, 4-20mA) and digital (RS-232) outputs, allowing for easy integration into automation systems, laboratory data acquisition setups, and industrial control platforms. These flexible output options enable real-time monitoring and remote data logging.

### Compact, Durable, and Reliable

Designed for tabletop use, the Model O2 is compact and lightweight (7.5  $\times$  3.0  $\times$  11.0 inches, 5 lbs) while maintaining industrial-grade durability. With no moving parts, it resists mechanical vibration and environmental disturbances, ensuring long-term reliability in demanding research and industrial environments.

# **Accessories**



#### Calibration Kit


Two Regulator Valves, two tanks of calibrating gases (21.00 and 99.99% +/- 0.05%) with Cal Kit Tubing Assembly and hard plastic carry case. (35 Liter bottles: approx 100 calibrations.

# Sensor Inlet Filter



PTFE moisture barrier/dust barrier for sensor, no fittings . (Package of 5). 25MM (package of 1)

Percent Oxygen     Percent CO2     Pressre Oxygen     Cell Temp Oxygen     Cell Temp Oxygen     Cell Temp Co2     Sample Flow     Marm Status			Calbrate High Calbrate High 100.00 Set 0 Calbrate Low 2010 Set 02 002 Cal High
Pressure Oxygen     Cell Temp Oxygen     Cell Temp Oxygen     Cell Temp Oxygen     Adam Status		Set Filters 6 O2 Filter 0 CO2 Filter Later Temp 33.84 °C	Calbrate High
Cell Temp Oxygen     Cell Temp Oxygen     Cell Temp CO2     Sample Flow     Marm Status		6 02 Filter 0 C02 Filter Laser Temp 33.84 °C	100.00 Set 0 Calibrate Low 20.00 Set 02
Cell Temp CO2     Sample Flow     Morm Status		6 02 Filter 0 C02 Filter Laser Temp 33.84 °C	100.00 Set 0 Calibrate Low 20.00 Set 02
Sample Flow     Alarm Status		0 CO2 Filter	Calibrate Low 20.00 Set 05
<ul> <li>Marm Status</li> </ul>		Laser Temp 33.84 °C	20.00 Set 02
		33.84 °C	
			C02 Cal High
- M 02			
	2	3.947 mA	0.00 Set CO
10 02 Pres (mb)		Cell Null Value	CO2 Cal Low
		298	32.50 Set CO
		Prediction Ga	
45.00 02 0 116		2 100 Set 0	02 Prediction Gain
Temp 0.00 002 POOP 13 21.00401.00401 Pred Threshold			
Log Data		5 100 Set C	02 Pred Threshold
J A C		C02 Pred Ga	m_l
		2 0 Set C	202 Prediction Gain
ot A Parto		a O Set C	
	0 Sample Row 0 Alarm 15.00 02 0 110 0.000 002 NORM 1210 Destroyed 1 A Petro	W Sample Row         0         Narm Status           45.00         02         0         1100         422           0.00         002         WORM 13.11.0446.0446         0442           Log Data           1         1         1	Organization         April 2000         April



#### OxiSoft Software

Oxisoft is a graphical oriented tool for controlling and displaying Oxigraf oxygen analyzers. Data logging of concentration, temperature and pressure for both oxygen and CO2. Dual panels. One panel for Data Display and one for setup and configuration.

# General Research Tabletop Oxygen Analyzer Item # 07-0006 Model O2

# **Technical Data**

Performance Conditions	
Ambient Temperature (Operating)	5 to 40 °C 40 to 102 °F
Ambient Temperature (Storage)	-20 to 60 °C -2 to 140 °F
Cell Pressure	10.2 to 17.4 psia 70 to 120 kPascal 500 to 900 mmHg
Warm-up for Full Accuracy	5 min
Altitude	Two point calibration required after change in altitude of 2000 feet
Humidity	0 to 95% non-condensing
Performance Specifications	
Range	5 to 100%
Modes	XC/LN Modes
Pump	Yes
Inlet Fittings	Luer
Resolution	0.1%
Accuracy - Stability (8 Hrs)	$\pm 0.2\%$ after 5 minute warm up in LN mode (nitrogen mixtures) $\pm 0.4\%$ after 5 minute warm up in XC mode
Flow (Using Pump)	50 to 250 ml/min adjustable.
Response Time	150 ms @ 250 ml/min flow. (12 software averaging filters available from front panel).
Analog Output	0 to 1.0 volts for 0 to 100% oxygen, 1.00K - 1% output resistance
Digital Output	RS232: 9600 baud default, 8 bit, no parity
Electrical Specifications	
Power Requirements	Voltage (DC)- 12 V Current- 1.5 A
External Power Supply	95 to 250 VAC, 47 to 63 Hz
Mechanical Specifications	
Dimensions (W x H x D)	7.5 x 3.0 x 11.0 in190 x 76 x 280 mm
Weight Instrument	Instrument 5 lbs (2.3Kg), power module 1.5 lbs (0.7Kg)

